

Human TSH (Thyroid Stimulating Hormone) Sandwich ELISA Kit Ready-To-Use (STJE0019225)

STJE0019225

GENERAL INFORMATION

Product Type	Kits
Short Description	This TSH Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	0.312-20mIU/L
Sensitivity	0.121mIU/L
Detection Limit	0.312-20mIU/L
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum.
Gene Symbol	
Uniprot ID	