

Mouse TSHR (Thyroid Stimulating Hormone Receptor) Sandwich ELISA Kit Ready-To-Use (STJE0019229)

STJE0019229

GENERAL INFORMATION

Product Type	Kits
Short Description	This TSHR Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection Sensitivity	62.5-4000pg/mL
Detection Limit	25.6pg/mL
Assay Time	62.5-4000pg/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates or other biological fluids.
Gene	Tshr
Symbol	
Uniprot ID	TSHR_MOUSE