

Human VDAC1 (Voltage Dependent Anion Channel Protein 1) Sandwich ELISA Kit Ready-To-Use (STJE0019380)

STJE0019380

GENERAL INFORMATION

Product Type	Kits
Short Description	This VDAC1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Sensitivity	31.2-2000pg/mL
Detection Limit	13.8pg/mL
Assay Time	31.2-2000pg/mL
Storage Instruction	3 hrs
Note	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates or other biological fluids.
Gene Symbol	VDAC1
Uniprot ID	VDAC1_HUMAN