

Human VMO1 (Vitelline Membrane Outer Layer Protein 1 Homolog (VMO1) ELISA Kit) Sandwich ELISA Kit Ready-To-Use (STJE0019448)

STJE0019448

GENERAL INFORMATION

Product Type	Kits
Short Description	This VMO1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection Sensitivity	78.1-5000pg/mL
Detection Limit	28.8pg/mL
Assay Time	78.1-5000pg/mL
Storage	3 hrs
Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Gene	VMO1
Symbol	
Uniprot ID	VMO1_HUMAN