

Mouse WNT5A (Wingless Type MMTV Integration Site Family, Member 5A) Sandwich ELISA Kit Ready-To-Use (STJE0019498)

STJE0019498

GENERAL INFORMATION

Product Type	Kits
Short Description	This WNT5A Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	31.2-2000pg/mL
Sensitivity	12.4pg/mL
Detection Limit	31.2-2000pg/mL
Assay Time	3 hrs
Storage Instruction	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene	Wnt5a
Symbol	
Uniprot ID	WNT5A_MOUSE