

2019-nCoV PL-Pro protein (His-Tag) (STJP000002)

STJP000002

GENERAL INFORMATION

Product Type	Proteins
Short Description	2019-nCoV PL-Pro-His-Tag protein was developed from e.coli and has a target region of His-Tag. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	Lyophilized
Formulation	
Purification	
Dilution Range	
Storage Instruction	Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081