

## 2019-nCoV S1 Protein protein (His-Tag) (STJP000005)

STJP000005

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | 2019-nCoV S1 Protein-His-Tag protein was developed from mammalian cells and has a target region of His-Tag. For use in research applications. |
| <b>Host/Source</b>       | Mammalian cells   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Concentration</b>       | Lyophilized  |
| <b>Formulation</b>         |  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Storage Instruction</b> | Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt. |

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Sequence**