

2019-nCoV Nucleocapsid protein protein (His-Tag) (STJP000006)

STJP000006

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | 2019-nCoV Nucleocapsid protein-His-Tag protein was developed from e.coli and has a target region of His-Tag. For use in research applications. |
| Host/Source | E.coli |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | Lyophilized |
| Formulation | |
| Purification | |
| Dilution Range | |
| Storage Instruction | Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt. |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence