

2019-nCoV S1 Protein Fc Tag protein (Fc-Tag) (STJP000007)

STJP000007

GENERAL INFORMATION

Product Type	Proteins
Short Description	2019-nCoV S1 Protein Fc Tag-Fc-Tag protein was developed from mammalian cells and has a target region of Fc-Tag. For use in research applications.
Host/Source	Mammalian cells

PRODUCT PROPERTIES

Concentration	Lyophilized
Formulation	
Purification	
Dilution Range	
Storage	Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence