

2019-nCoV Spike protein fragment 2 protein (His-Tag) (STJP000009)
STJP000009**GENERAL INFORMATION**

Product Type	Proteins
Short Description	2019-nCoV Spike protein fragment 2-His-Tag protein was developed from e.coli and has a target region of His-Tag. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	Lyophilized
Formulation	
Purification	
Dilution Range	
Storage	Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence