

## 2019-nCoV 3C-like Proteinase protein (GST-Tag) (STJP000012)

STJP000012

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	2019-nCoV 3C-like Proteinase-GST-Tag protein was developed from e.coli and has a target region of GST-Tag. For use in research applications.
<b>Host/Source</b>	E.coli

### PRODUCT PROPERTIES

<b>Concentration</b>	Lyophilized
<b>Formulation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Storage Instruction</b>	Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt.

### TARGET INFORMATION

Gene ID  
Gene Symbol  
Uniprot ID  
Immunogen  
Sequence