

2019-nCoV NSP9 protein (His-Tag) (STJP000016)

STJP000016

GENERAL INFORMATION

Product Type	Proteins
Short Description	2019-nCoV NSP9-His-Tag protein was developed from e.coli and has a target region of His-Tag. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	Lyophilized
Formulation	
Purification	
Dilution Range	
Storage Instruction	Store at 2-8 °C for one week. Thereafter store at -20 to -80 °C for 1 year from the date of receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	