

APN Protein (Mammalian, C-6His) (STJP000049)
STJP000049**GENERAL INFORMATION**

Product Type	Proteins
Short Description	APN Protein (Mammalian, C-6His) (STJP000049) produced in human.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 mg/ml filtered solution of PBS, pH7.4.
Purification	Greater than 95% as determined by reducing SDS-PAGE
Dilution Range	
Storage Instruction	Store at -20°C. Stable for 6 months. Minimize freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen Sequence	