

Recombinant 2019-nCoV 3C-like Proteinase Protein (His-Tag) (Ser1-Gln306) (STJP000060)

STJP000060

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV 3C-like Proteinase Protein His-Tag is developed from E.coli and has a target region of Ser1-Gln306. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of PBS, 1 mM EDTA, 10% Glycerol, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	