

Recombinant 2019-nCoV NSP1 Protein (His-Tag) (Met1-Gly180) (STJP000062)

STJP000062

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV NSP1 Protein His-Tag is developed from E.coli and has a target region of Met1-Gly180. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of 20 mM Tris-HCl, 150 mM NaCl, 10% Glycerol, pH 8.5.
Purification	
Dilution Range	> 85 % as determined by reducing SDS-PAGE.
Storage	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	