

## Recombinant 2019-nCoV Papain-Like Protease Protein (Glu1564-Lys1878) (STJP000067)

STJP000067

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV Papain-Like Protease Protein NA is developed from E.coli and has a target region of Glu1564-Lys1878. For use in research applications.
<b>Host/Source</b>	E.coli

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Supplied as a 0.2 $\mu$ M filtered solution of 20mM Tris-HCl, 10 mM Mercaptoethanol [BME], 20% Glycerol, pH 7.5.
<b>Purification</b>	
<b>Dilution Range</b>	> 90 % as determined by reducing SDS-PAGE.
<b>Storage</b>	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Sequence</b>	