

Recombinant 2019-nCoV Spike Protein (RBD-SD1, mFc-Tag) (Arg319-Pro589) (STJP000071)

STJP000071

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV Spike Protein RBD-SD1, mFc-Tag is developed from Human and has a target region of Arg319-Pro589.
Description	For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	