

## Recombinant 2019-nCoV Spike Protein (RBD-SD1, Avi-His-Tag) (Arg319-Ser591) {Biotin} (STJP000073)

STJP000073

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV Spike Protein RBD-SD1, Avi-His-Tag is developed from Human and has a target region of Arg319-Ser591. For use in research applications.
<b>Host/Source</b>	Human

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Storage Instruction</b>	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Sequence</b>	