

## Recombinant 2019-nCoV S1 Protein (Fc-Tag) (Gln14-Arg685) (STJP000075)

STJP000075

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV S1 Protein Fc-Tag is developed from Human and has a target region of Gln14-Arg685. For use in research applications.
<b>Host/Source</b>	Human

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	> 90 % as determined by reducing SDS-PAGE.
<b>Storage Instruction</b>	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Sequence**