

Recombinant 2019-nCoV S1 Protein (His-Tag) (Val16-Arg685) (STJP000077)

STJP000077

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV S1 Protein His-Tag is developed from Human and has a target region of Val16-Arg685. For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage Instruction	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence