

Recombinant 2019-nCoV S1 Protein (CTD, His-Tag) (STJP000078)

STJP000078

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | Recombinant 2019-nCoV S1 Protein CTD, His-Tag is developed from Human. For use in research applications. |
| Host/Source | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | Supplied as a 0.2 μ M filtered solution of PBS, pH 7.4. |
| Purification | |
| Dilution Range | > 95 % as determined by reducing SDS-PAGE. |
| Storage Instruction | Store at $<-20^{\circ}\text{C}$, stable for 6 months. Please minimize freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen Sequence | |