

## Recombinant 2019-nCoV Spike Protein (NTD, His-Tag) (STJP000079)

STJP000079

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV Spike Protein NTD, His-Tag is developed from Human. For use in research applications.
<b>Host/Source</b>	Human

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Supplied as a 0.2 $\mu$ M filtered solution of PBS, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Storage Instruction</b>	Store at $-20^{\circ}\text{C}$ , stable for 6 months. Please minimize freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen Sequence</b>	