

Recombinant 2019-nCoV NP NTD domain (His-Tag) (STJP000081)

STJP000081

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV NP NTD domain His-Tag is developed from E.coli. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of PBS, 2M Urea, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage Instruction	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen Sequence	