

Recombinant 2019-nCoV Envelope Protein (His-Tag) (Met1-Val75) (STJP000083)

STJP000083

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV Envelope Protein His-Tag is developed from E.coli and has a target region of Met1-Val75. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of 20 mM Tris-HCl, 200 mM NaCl, pH 8.0.
Purification	
Dilution Range	> 85 % as determined by reducing SDS-PAGE.
Storage	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence