

Recombinant Human ACE2 Protein (Fc-Tag) (Gln18-Ser740) (STJP000086)

STJP000086

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant Human ACE2 Protein Fc-Tag is developed from Human and has a target region of Gln18-Ser740. For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of 20 mM Tris-HCl, 300 mM NaCl, 1 mM ZnCl ₂ , 10% Glycerol, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage Instruction	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence