

## Recombinant Human ACE2 Protein (Avi-His-Tag) (Gln18-Ser740) (STJP000087)

STJP000087

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant Human ACE2 Protein Avi-His-Tag is developed from Human and has a target region of Gln18-Ser740. For use in research applications.
<b>Host/Source</b>	Human

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Supplied as a 0.2 µM filtered solution of 20 mM Tris-HCl, 300 mM NaCl, 1 mM ZnCl <sub>2</sub> , 10% Glycerol, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Storage</b>	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Sequence</b>