

Recombinant Human APN Protein (His-Tag) (Lys69Lys967) (STJP000089)

STJP000089

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant Human APN Protein His-Tag is developed from Human and has a target region of Lys69Lys967. For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage Instruction	Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence