

## Recombinant Mouse APN Protein (His-Tag) (Lys69Ser966) (STJP000090)

STJP000090

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | Recombinant Mouse APN Protein His-Tag is developed from Human and has a target region of Lys69Ser966. For use in research applications. |
| <b>Host/Source</b>       | Human   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Concentration</b>       |   |
| <b>Formulation</b>         | Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.                    |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | > 95 % as determined by reducing SDS-PAGE.                                |
| <b>Storage Instruction</b> | Store at <-20°C, stable for 6 months. Please minimize freeze-thaw cycles. |

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Sequence**