

Recombinant 2019-nCoV Spike Protein (RBD, His-Tag) (Arg319-Phe541) {Biotin} (STJP000095)

STJP000095

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV Spike Protein RBD, His-Tag is developed from Baculovirus-Insect Cells and has a target region of Arg319-Phe541. For use in research applications.
Host/Source	Baculovirus-Insect Cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from sterile PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8 °C for 2-7 days. Aliquots of reconstituted samples are stable at <-20 °C for 3 months.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	