

## Recombinant 2019-nCoV S2 Protein (ECD, His-Tag) (Ser686-Pro1213) (STJP000099)

STJP000099

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV S2 Protein ECD, His-Tag is developed from Baculovirus-Insect Cells and has a target region of Ser686-Pro1213. For use in research applications.
<b>Host/Source</b>	Baculovirus-Insect Cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilized from sterile 20 mM PB, 300 mM NaCl, 10 % glycerol, pH 7.0.
<b>Purification</b>	
<b>Dilution Range</b>	> 90 % as determined by reducing SDS-PAGE.
<b>Storage</b>	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8 °C
<b>Instruction</b>	for 2-7 days. Aliquots of reconstituted samples are stable at <-20 °C for 3 months.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Sequence</b>	