

Recombinant 2019-nCoV Spike Protein (RBD, Fc-Tag) (Arg319-Phe541) (STJP000112)

STJP000112

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV Spike Protein RBD, Fc-Tag is developed from HEK293 and has a target region of Arg319-Phe541. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from sterile PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE. > 95 % as determined by SEC-HPLC.
Storage	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8 °C
Instruction	for 2-7 days. Aliquots of reconstituted samples are stable at <-20 °C for 3 months.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	