

Recombinant 2019-nCoV Spike Protein (RBD, rFc-Tag) (Arg319-Phe541 & V367F) (STJP000117)

STJP000117

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV Spike Protein RBD, rFc-Tag is developed from HEK293 and has a target region of Arg319-Phe541 & V367F.
Description	For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from sterile PBS, pH 7.4.
Purification	
Dilution Range	> 90 % as determined by reducing SDS-PAGE.
Storage	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8 °C
Instruction	for 2-7 days. Aliquots of reconstituted samples are stable at <-20 °C for 3 months.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence