

## Recombinant 2019-nCoV Spike Protein (RBD, His-Tag) (Arg319-Phe541 & A435S) (STJP000122)

STJP000122

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV Spike Protein RBD, His-Tag is developed from HEK293 and has a target region of Arg319-Phe541 & A435S.
<b>Description</b>	For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilized from sterile PBS, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Storage</b>	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8 °C
<b>Instruction</b>	for 2-7 days. Aliquots of reconstituted samples are stable at <-20 °C for 3 months.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Sequence</b>	