

Recombinant 2019-nCoV S1 Protein (RBD, His-Tag) (Val16-Arg685 & D614G) (STJP000124)

STJP000124

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV S1 Protein RBD, His-Tag is developed from HEK293 and has a target region of Val16-Arg685 & D614G.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from sterile PBS, pH 7.4
Purification	
Dilution Range	> 90 % as determined by reducing SDS-PAGE.
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	