

## Recombinant 2019-nCoV S1 Protein (Avi-His-Tag) (Val16-Arg685) {Biotin} (STJP000125)

STJP000125

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant 2019-nCoV S1 Protein Avi-His-Tag is developed from HEK293 and has a target region of Val16-Arg685. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilized from sterile PBS, pH 7.4
<b>Purification</b>	
<b>Dilution Range</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Storage</b>	
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Sequence</b>	