

Recombinant 2019-nCoV Spike Protein (RBD, His-Tag) (Arg319-Phe541 & W436R) (STJP000127)

STJP000127

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV Spike Protein RBD, His-Tag is developed from HEK293 and has a target region of Arg319-Phe541 & W436R. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from sterile PBS, pH 7.4
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	