

## Recombinant 2019-nCoV Spike Protein (RBD, His-Tag) (Arg319-Phe541) {Biotin} (STJP000129)

STJP000129

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | Recombinant 2019-nCoV Spike Protein RBD, His-Tag is developed from HEK293 and has a target region of Arg319-Phe541. For use in research applications. |
| <b>Host/Source</b>       | HEK293  |

### PRODUCT PROPERTIES

|                       |  |
|-----------------------|--|
| <b>Concentration</b>  |  |
| <b>Formulation</b>    | Lyophilized from sterile PBS, pH 7.4       |
| <b>Purification</b>   |  |
| <b>Dilution Range</b> | > 95 % as determined by reducing SDS-PAGE. |
| <b>Storage</b>        |  |
| <b>Instruction</b>    |  |

### TARGET INFORMATION

|                    |  |
|--------------------|--|
| <b>Gene ID</b>     |  |
| <b>Gene Symbol</b> |  |
| <b>Uniprot ID</b>  |  |
| <b>Immunogen</b>   |  |
| <b>Sequence</b>    |  |