

Recombinant 2019-nCoV NP (His-Tag) (STJP000130)

STJP000130

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV NP His-Tag is developed from E.coli. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 μ M filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage Instruction	Store at $<-20^{\circ}\text{C}$, stable for 7 months. Please minimize freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen Sequence	