

Recombinant 2019-nCoV NP NTD domain V2 (His-Tag) (Gly44-Ser180) (STJP000131)

STJP000131

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV NP NTD domain V2 His-Tag is developed from E.coli and has a target region of Gly44-Ser180. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by reducing SDS-PAGE.
Storage	Store at <-20°C, stable for 8 months. Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	