

Recombinant HCoV-OC43 S1+S2 Protein (ECD, His-Tag) (Met1-Pro1304) (STJP000132)

STJP000132

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant HCoV-OC43 S1+S2 Protein ECD, His-Tag is developed from Baculovirus-Insect Cells and has a target region of Met1-Pro1304. For use in research applications.
Host/Source	Baculovirus-Insect Cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from sterile 20 mM PB, 300 mM NaCl, 10 % glycerol, pH 7.5
Purification	
Dilution Range	> 85 % as determined by reducing SDS-PAGE.
Storage	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8 °C
Instruction	for 2-7 days. Aliquots of reconstituted samples are stable at <-20 °C for 3 months.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	