

Recombinant SARS-CoV-2 S protein RBD (C-mFc-6His-tag) (Arg319-Phe541) (STJP000135)

STJP000135

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant SARS-CoV-2 S protein RBD C-mFc-6His-tag is developed from HEK293 and has a target region of Arg319-Phe541.
Description	For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95 % as determined by SDS-PAGE.
Storage	Store at <-20°C, stable for 6 months after receipt. NA NA Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	