

Recombinant 2019-nCoV NP CTD domain V2 (C-6His) (Thr247-Pro364) (STJP000137)

STJP000137

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant 2019-nCoV NP CTD domain V2 C-6His is developed from E.coli and has a target region of Thr247-Pro364. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.2 µM filtered solution of 50mM PB, 150mM NaCl, 1mM EDTA, pH 7.4
Purification	
Dilution Range	>95% as determined by reducing SDS-PAGE.
Storage	Store at <-20°C, stable for 6 months after receipt. NA NA Please minimize freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	