

CDNF protein (STJP000138)

STJP000138

GENERAL INFORMATION

Product Type Proteins
Short Description CDNF-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation PBS pH 7.4
Purification Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media. Protein is sterile-filtrated through 0.22 µm filter
Dilution Range >90%
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [441549](#)
Gene Symbol [CDNF](#)
Uniprot ID [CDNF_HUMAN](#)
Immunogen
Sequence