

MANF protein (STJP000139)

STJP000139

GENERAL INFORMATION

Product Type Proteins
Short Description MANF-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation PBS pH 7.4, concentration
Purification Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID 7873
Gene Symbol MANF
Uniprot ID MANF_HUMAN
Immunogen Sequence