

## 2B4 protein (STJP000141)

STJP000141

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	2B4-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	PBS pH 6.0
<b>Purification</b>	Purified by Ni-affinity chromatography, ion-exchange chromatography and desalting chromatography from serum-free CHO growth media
<b>Dilution Range</b>	
<b>Storage Instruction</b>	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

<b>Gene ID</b>	18106
<b>Gene Symbol</b>	Cd244
<b>Uniprot ID</b>	CD244_MOUSE
<b>Immunogen</b>	
<b>Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081