

CD48 protein (STJP000142)

STJP000142

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD48-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Purified by ion-exchange chromatography, hydrophobic interaction chromatography and gel-filtration from serum-free CHO growth media
Dilution Range	
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	12506
Gene Symbol	Cd48
Uniprot ID	CD48_MOUSE
Immunogen	
Sequence	