

BDNF protein (STJP000143)

STJP000143

GENERAL INFORMATION

Product Type Proteins
Short Description BDNF-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation PBS pH 7.4
Purification Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media. Protein is sterile-filtrated through 0.22 µm filter
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [627](#)
Gene Symbol [BDNF](#)
Uniprot ID [BDNF_HUMAN](#)
Immunogen
Sequence