

## Recombinant Bovine DNase I protein (STJP000145)

STJP000145

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** Recombinant Bovine DNase I-protein was developed for use in research applications.  
**Host/Source**

### PRODUCT PROPERTIES

**Concentration** 0.1 mg/ml  
**Formulation** 10mM Ca-acetate pH 5.0, 140 mM NaCl, 5mM MgCl<sub>2</sub>  
**Purification** Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media  
**Dilution Range**  
**Storage Instruction** Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen Sequence**