

## DNase I protein (STJP000146)

STJP000146

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** DNase I-protein was developed for use in research applications.  
**Host/Source**

### PRODUCT PROPERTIES

**Concentration** 0.5 mg/ml  
**Formulation** PBS pH 7.4  
**Purification** Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media  
**Dilution Range**  
**Storage Instruction** Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

**Gene ID** [1773](#)  
**Gene Symbol** [DNASE1](#)  
**Uniprot ID** [DNASE1\\_HUMAN](#)  
**Immunogen Sequence**