

**GDNF protein (STJP000148)**  
STJP000148**GENERAL INFORMATION**

**Product Type** Proteins  
**Short Description** GDNF-protein was developed for use in research applications.  
**Host/Source**

**PRODUCT PROPERTIES**

**Concentration** 0.96 mg/ml  
**Formulation** PBS pH 7.4  
**Purification** Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media. Protein is sterile-filtrated through 0.22 µm filter  
**Dilution Range**  
**Storage Instruction** Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

**TARGET INFORMATION**

**Gene ID** [2668](#)  
**Gene Symbol** [GDNF](#)  
**Uniprot ID** [GDNF\\_HUMAN](#)  
**Immunogen**  
**Sequence**